

Sprint S237

MallMASTER

Compact/portable scooter

The Eclipse Mallmaster series of compact scooters offer sleek designs and detachable batteries for ultimate lightweight portability. This class of scooters are ideal for individuals going shorter distances and needing to transport their unit to get from place to place.

- > Mallmaster Sprint 3-wheel scooter
- > Lightweight and portable
- > 22AH detachable battery pack with 1.8A charger
- > Flat free tires
- > Deluxe swivel seat with sliders
- > Colours: red, blue or white



\$2,295.00
MSRP



S237 Specifications

Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	Driving dist.	Turning radius	Safe climb°
136KG/300LBS	40" x 21" x 32"	47KG/104LBS	0.6HP	7.2KPH/4.5MPH	2 x 22AH	F8" R9"	17KM/10.6M	86CM/34"	6°/10%

All specifications subject to change without notice.

Dealer stamp



Discovering your world again—in style

www.eclipsemedical.com

#112019EN