

# Trailblazer 889<sup>SL/SE</sup> SHOPPRIDER

## Luxury scooter

Shoprider's Crown Trailblazer 889<sup>SL/SE</sup> provides feature rich luxury in personal mobility. Experience the outdoors in unprecedented style, comfort and performance. The Crown Series are the most popular scooters in the Shoprider Canada portfolio.

- > **Special Edition** Shoprider 4-wheel scooter
- > Delta tiller with 2-speed (hi-lo)
- > 2HP motor with 13.5 inch tires
- > Full suspension, full light package
- > 50AH batteries with 4A charger
- > Available colours: red, black or silver



### 889<sup>SL/SE</sup>

Specs	Weight capacity	Dimensions L x W x H	Total wt. w/bat.	Motor	Speed	Battery	Tire	24v charger	Driving dist.	Turning radius	Safe climb <sup>o</sup>
889 <sup>SL/SE</sup>	158KG/350LBS	49" x 24" x 45"	100KG/220LBS	2.0HP	15KPH/9.45MPH	2 x 50AH	4 x 13.5"	4AMP/HR	40KM/25M	1500MM/59"	9°/15%

All specifications subject to change without notice.



**MOOSE MOBILITY SCOOTER CORP.**  
 7113 Argyll Road Edmonton AB T6C 4A5  
 bus: 780.465.0856  
 toll free: 1.844.MYMOOSE  
 info@moosemobilityscooter.com



Discovering your world again—in style

[www.eclipsemedical.com](http://www.eclipsemedical.com)

#092013EN