

Item #	Product Dimensions	UOM
Maverick20CS	65" x 31" x 56	1/ea
Maverick22CS	65" x 31" x 56	1/ea



Maverick Executive Scooter

- Sport Style Captain's Seat available in 20" or 22" seat widths.
- Captain's seat swivels and features a reclining backrest, and height adjustable headrest.
- Padded armrests feature easy turn cylinder for angle adjustment.
- 19.5" Air-filled, low-profile tires. (Figure E)
- Front and rear motorcycle inspired suspension for the smoothest outdoor ride quality. (Figure G)
- Stylish controls are conveniently located on the handlebars for ease of use. (Figure F)
- Smoked windshield provides glare reduction on the Digital LED display with real-time feedback. (Figure D)
- Efficient, extra bright LED lighting. (Figure C)
- Easy to access key switch and charge port located at the control panel. (Figure F)
- Hazard lights, turn signals and reflectors for superb visibility.
- Rear view mirror included.
- Positioning belt included.
- Anti-tip wheels.
- Easy, freewheel mode.

specifications

Weight Capacity	400 lbs.	Charger	8A
Top Speed	9.5 mph	Brakes	Electromagnetic & Hand Brake
Max Range	35 miles	Front Wheels	19.5" x 2.75"
Turning Radius	62"	Rear Wheels	19.5" x 2.75"
Climbing Angle	10°	Base Weight	181 lbs.
Ground Clearance	5"	Battery Weight	94 lbs.
Dimensions	65" x 31"	Seat Weight	55 lbs.
Seat Dimensions	20" x 20" ; 22" x 20"	Total Weight	330 lbs.
Controller	Dynamic 160A	Warranty on Frame	Limited Lifetime
Motor	1300W	Warranty on Electronic Controller/Drive Train Components	14 months
Batteries	75 AH	Warranty on Battery	6 months