

POWERCHAIR**SPECIFICATIONS**

<i>Model</i>	CTM HS-1500	CTM HS-2850	SupaChair Safari
<i>Weight Capacity</i>	115.67kg/255lbs	136kg/300lbs	150kg/330lbs
<i>Dimensions (LxWxH)</i>	31.9"x23.6"x36.8"	39"x24"x42.5"	34.4"x23.6"x33"
<i>Seat Width</i>	18.5"	-	-
<i>Ground Clearance</i>	1.6"	2.4"	-
<i>Total Weight with Battery</i>	47.81kg/105.4lbs	82.1kg/181lbs	36.5kg/80lbs
<i>Speed</i>	6.44kph/4mph	6.76kph/4.2mph	6.2kph/3.9mph
<i>Drive Wheel</i>	8" (rear) 5" (castor)	10" (middle) 6" (castor)	8" (rear)
<i>Battery</i>	2 x 12V 14 ah	2 x 12V 36 ah	24V 15 ah (SLA)
<i>Charger</i>	2 amp offboard	5 amp offboard	2 amp offboard
<i>Motor</i>	150W	420W	-
<i>Driving Distance</i>	8.37km/5.2miles	19.96km/12.4miles	10km/6.2miles
<i>Turning Radius</i>	695.96mm/27.4"	619.76mm/24.4"	700mm/27.6"
<i>Safe Climb°</i>	6°	6°	9°
<i>Suspension</i>	-	Front, rear, and centre	Dual seat post

"-" indicates that the information is unavailable at this time. Specifications are subject to change without notice.